

MMP-8 / MMP8 Antibody (Biotin), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10254-R102-B

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human MMP-8 / MMP8 Protein (Catalog#10254-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human MMP-8 / MMP8 (rh MMP-8 / MMP8; Catalog#10254-H08H; NP_002415.1; Met1-Gly467), and then biotinylated.
Ig Type:	Rabbit IgG
Clone ID:	102
Specificity:	Human MMP-8 / MMP8
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/ml

Please Note: Optimal concentrations/dilutions should be determined by the end user.