

CD64/FCGR1 Antibody (PerCP), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10256-R401-C

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human CD64/FCGR1 Protein (Catalog#10256-H08H)
Reagents:	PerCP-conjugated Rabbit monoclonal antibody
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD64/FCGR1 (rh CD64/FCGR1; Catalog#10256-H08H; NP_000557.1; Met1-Pro288) and conjugated with PerCP under optimum conditions, the unreacted PerCP was removed.
Ig Type:	Rabbit IgG
Clone ID:	401
Specificity:	Human CD64/FCGR1
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2℃-8℃ for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !
Alternative Names:	FCGR1

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.


CD64/FCGR1 Antibody (PerCP), Rabbit MAb

Catalog Number: 10256-R401-C



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Flow cytometric analysis of Human CD64 expression on human whole blood monocytes. Cells were stained with PerCP-conjugated anti-Human CD64. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable monocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.