

AXL Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10279-MM17T-P

GENERAL INFORMATION

Immunogen:	Recombinant Human AXL Protein (Catalog#10279-H08H)
Reagents:	PE-conjugated Mouse MAb antibody
Preparation	This product is a recombinant monoclonal antibody expressed from HEK293 cells and conjugated with PE under optimum conditions, the unreacted PE was removed.
Ig Type:	Mouse IgG1
Clone ID:	17
Specificity:	Human AXL
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

AXL Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10279-MM17T-P



Flow cytometric analysis of Human AXL expression on DU145 cells. Cells were stained with PE-conjugated anti-Human AXL. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>