

CCSAP Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 102870-T08

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the N-terminus of the Human CCSAP
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human CCSAP, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CCSAP
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	CCSAP

APPLICATIONS

Applications:	IHC-P
---------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:100-1:500
-------	--------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CCSAP Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 102870-T08



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human CCSAP in human brain with rabbit polyclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.