

**GDI1 Antibody, Rabbit PAb,
Antigen Affinity Purified**



Sino Biological
Biological Solution Specialist

Catalog Number: 103113-T32

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the N-terminus of the Human GDI1
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human GDI1, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Mouse,RAT,MACFA, (Species predicted to react based on 100% sequence homology)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB
----------------------	----

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:2000
---------------------	------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

GD11 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 103113-T32



Anti-GD11 rabbit polyclonal antibody at 1:500 dilution
Lane A: Hela Whole Cell Lysate
Lane B: SH-SY5Y Whole Cell Lysate

Lysates/proteins at 30 µg per lane.
Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 51 kDa
Observed band size: 51 kDa