

Fibronectin Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10314-T16

GENERAL INFORMATION

Immunogen:	10314-H08H
Preparation	Produced in rabbits immunized with purified, recombinant Human Fibronectin (rh Fibronectin; Catalog#10314-H08H; NP_997647.1; Ser607-Pro1265). Fibronectin specific IgG was purified by Human Fibronectin affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Fibronectin
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	CIG

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
-------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.