

HGF Activator/HGFA Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10329-T16

GENERAL INFORMATION

Immunogen:	Recombinant Human HGF Activator/HGFA Protein (Catalog#10329-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human HGF Activator/HGFA (rh HGF Activator/HGFA; Catalog#10329-H08H; NP_001519.1; Met1-Ser655). HGF Activator/HGFA specific IgG was purified by Human HGF Activator/HGFA affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human HGF Activator/HGFA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
--------------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.