

# ICAM-2/CD102 Antibody, Rabbit MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 10332-R022

EliteRmab® is a registered trademark of Sino Biological Inc.

| GENERAL INFORMATION            |                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen:                     | Recombinant Human ICAM-2/CD102 Protein (Catalog#10332-H08H)                                                                                                                                                                                                                                                                                                                                         |
| Preparation                    | This antibody was obtained from a rabbit immunized with purified, recombinant Human ICAM-2/CD102 (rh ICAM-2/CD102; Catalog#10332-H08H; NP_000864.2; Met1-Gln223).                                                                                                                                                                                                                                   |
| Ig Type:                       | Rabbit IgG                                                                                                                                                                                                                                                                                                                                                                                          |
| Clone ID:                      | 022                                                                                                                                                                                                                                                                                                                                                                                                 |
| Specificity:                   | Human ICAM-2/CD102                                                                                                                                                                                                                                                                                                                                                                                  |
| Formulation:                   | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                                                                                                                                                                     |
| Storage:                       | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.<br>Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles. |
| APPLICATIONS                   |                                                                                                                                                                                                                                                                                                                                                                                                     |
| Applications:                  | ELISA(Cap)                                                                                                                                                                                                                                                                                                                                                                                          |
| RECOMMENDED CONCENTRATION      |                                                                                                                                                                                                                                                                                                                                                                                                     |
| Sandwich ELISA<br>(Capture Ab) | ELISA(Cap): 1:250-1:2000<br>This antibody will detect Human ICAM-2 / CD102 in ELISA pair set (Catalog: # SEK10332). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 10332-R015).                                                                                                                                                                               |

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**