

ICAM-2 / CD102 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10332-RP02

GENERAL INFORMATION	
Immunogen:	10332-H08H
Preparation	Produced in rabbits immunized with purified, recombinant Human ICAM-2 / CD102 (rh ICAM-2 / CD102; Catalog#10332-H08H; NP_000864.2; Met 1-Gln 223). ICAM-2 / CD102 specific IgG was purified by Human ICAM-2 / CD102 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human ICAM-2 / CD102
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	CD102
APPLICATIONS	
Applications:	WB,ELISA,IP
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:500-1:1000
Immunoprecipitation	IP: 1-2 µL/mg of lysate
ELISA	ELISA: 1:25000-1:50000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

ICAM-2 / CD102 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10332-RP02



Anti-ICAM2 rabbit polyclonal antibody at 1:500 dilution

Lane A: HL60 Membrane Lysate

Lysates/proteins at 20 µg per lane.

Secondary


Goat Anti-Rabbit IgG H&L (Dylight800) at 1/10000 dilution.

Developed using the Odyssey technique.

Performed under reducing conditions.

Predicted band size: 31 kDa

Observed band size: 47 kDa



ICAM2 was immunoprecipitated using:

Lane A: 0.5 mg Hela Whole Cell Lysate

0.5 µL anti-ICAM2 rabbit polyclonal antibody and 60 µg of Immunomagnetic beads Protein G.

Primary antibody:

Anti-ICAM2 rabbit polyclonal antibody, at 1:500 dilution

Secondary antibody:

Dylight 800-labeled antibody to rabbit IgG (H+L), at 1:5000 dilution

Developed using the odssey technique.

Performed under reducing conditions.

Predicted band size: 31 kDa

Observed band size: 31 kDa