

IFNGR1 / CD119 / IFN-gamma-R1 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10338-R010

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD119 protein (Catalog#10338-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human IFNGR1 / CD119 (rh IFNGR1; Catalog#10338-H08H; NP_000407.1; Met 1-Gly 245).
Ig Type:	Rabbit IgG
Clone ID:	10
Specificity:	Human IFNGR1 / CD119
	No cross-reactivity in ELISA with Human IL10RB / CRFB4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human IFNGR1.

Please Note: Optimal concentrations/dilutions should be determined by the end user.