

Granzyme B/GZMB Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10345-T48

GENERAL INFORMATION	
Immunogen:	Recombinant Human Granzyme B/GZMB Protein (Catalog#10345-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human Granzyme B/GZMB (rh Granzyme B/GZMB; Catalog#10345-H08H; NP_004122.1; Met1-Tyr247). Granzyme B/GZMB specific IgG was purified by Human Granzyme B/GZMB affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Granzyme B/GZMB
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB, ELISA
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human Granzyme B/GZMB.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Granzyme B/GZMB Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10345-T48



Anti-GranzymeB rabbit polyclonal antibody at 1:500 dilution
Lane A: HL60 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.
Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 28 kDa
Observed band size: 28 kDa
(We are unsure as to the identity of these extra bands.)