

ECM1 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10362-MM03

GENERAL INFORMATION	
Immunogen:	Recombinant Human ECM1 protein (Catalog#10362-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human ECM1 (rhECM1; Catalog#10362-H08H; Met 1-Glu 540; NP_004416.1). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	6C3C6
Specificity:	Human ECM1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human ECM1.

Please Note: Optimal concentrations/dilutions should be determined by the end user.