

# IL-33 Antibody, Rabbit MAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 10368-R839

EliteRmab® is a registered trademark of Sino Biological Inc.

## GENERAL INFORMATION

|                     |                                                                                                                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen:</b>   | Recombinant Human IL-33 Protein (Catalog#10368-HNAE)                                                                                                                                                                                               |
| <b>Preparation</b>  | This antibody was obtained from a rabbit immunized with purified, recombinant Human IL-33 (rh IL-33; Catalog#10368-HNAE; NP_254274.1; Ser112-Thr270).                                                                                              |
| <b>Ig Type:</b>     | Rabbit IgG                                                                                                                                                                                                                                         |
| <b>Clone ID:</b>    | 839                                                                                                                                                                                                                                                |
| <b>Specificity:</b> | Human IL-33                                                                                                                                                                                                                                        |
| <b>Formulation:</b> | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                    |
| <b>Storage:</b>     | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles. |

## APPLICATIONS

|                      |              |
|----------------------|--------------|
| <b>Applications:</b> | ELISA, IHC-P |
|----------------------|--------------|

## RECOMMENDED CONCENTRATION


|              |                       |
|--------------|-----------------------|
| <b>IHC-P</b> | IHC-P: 1:500-1:2000   |
| <b>ELISA</b> | ELISA: 1:5000-1:10000 |

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**




Catalog Number: 10368-R839

EliteRmab® is a registered trademark of Sino Biological Inc.



Immunochemical staining of human IL-33 in human tonsil with rabbit monoclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.



Immunochemical staining of human IL-33 in human stomach with rabbit monoclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.