

CD32a/FCGR2A/Fc gamma RIIA Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10374-MM02-P

GENERAL INFORMATION

Immunogen:	Recombinant Human CD32a/FCGR2A/Fc gamma RIIA Protein (Catalog#10374-H08H)
Reagents:	PE-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD32a/FCGR2A/Fc gamma RIIA (rh CD32a/FCGR2A/Fc gamma RIIA; Catalog#10374-H08H; P12318-1; Met1-Ile218,167R) and conjugated with PE under optimum conditions, the unreacted PE was removed.
Ig Type:	Mouse IgG1
Clone ID:	02
Specificity:	Human CD32a/FCGR2A/Fc gamma RIIA
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !
Alternative Names:	FCGR2A

APPLICATIONS

Applications:	FCM
---------------	-----

RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD32a/FCGR2A/Fc gamma RIIA Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10374-MM02-P

Flow cytometric analysis of Human CD32a expression on human whole blood lymphocytes (left panel), monocytes (middle panel) and granulocytes (right panel). Cells were stained with PE-conjugated anti-Human CD32a. The histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes, monocytes and granulocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.