

CD16/Fc gamma RIII Antibody (PerCP), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10389-MM23-C

GENERAL INFORMATION

Immunogen:	Recombinant Human CD16a/Fc gamma RIIIa Protein (Catalog#10389-H08H)
Reagents:	PerCP-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD16a/Fc gamma RIIIa (rh CD16a/Fc gamma RIIIa; Catalog#10389-H08H; P08637-1; Met1-Gln208,176F) and conjugated with PerCP under optimum conditions, the unreacted PerCP was removed.
Ig Type:	Mouse IgG1
Clone ID:	23
Specificity:	Human CD16 (CD16a & CD16b)
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.09% sodium azide
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.
Alternative Names:	FCGR3

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD16/Fc gamma RIII Antibody (PerCP), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10389-MM23-C



Flow cytometric analysis of Human CD16 expression on human whole blood lymphocytes. Cells were stained with PerCP-conjugated anti-Human CD16. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of viable lymphocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.