

IL20RA / IL-20RA Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10397-R038

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human IL20RA / IL-20RA protein (Catalog#10397-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human IL20RA / IL-20RA (rh IL20RA / IL-20RA; Catalog#10397-H08H; NP_055247.3; Met 1-Lys 250).
Ig Type:	Rabbit IgG
Clone ID:	038
Specificity:	Human IL20RA / IL-20RA
	No cross-reactivity in ELISA with Human IFNGR1 Human IL10RB
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA, IHC-P
----------------------	--------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:100-1:500
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human IL20Ra.

Please Note: Optimal concentrations/dilutions should be determined by the end user.


IL20RA / IL-20RA Antibody, Rabbit MAb

Catalog Number: 10397-R038



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Immunohistochemical staining of human IL20Ra in human skin with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections). Positive staining was localized to stratified squamous epithelium.