

CD4 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10400-T24

GENERAL INFORMATION

Immunogen:	Recombinant Human CD4 Protein (Catalog#10400-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human CD4 (rh CD4; Catalog#10400-H08H; NP_000607.1; Met1-Trp390). CD4 specific IgG was purified by Human CD4 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CD4
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA, IHC-P
---------------	--------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD4 Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 10400-T24



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human CD4 in human tonsil with rabbit polyclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CD4 in human spleen with rabbit polyclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.