

Cystatin E / CST6 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10438-R022

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human Cystatin E / CST6 Protein (Catalog#10438-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human Cystatin E / CST6 (rh Cystatin E / CST6; Catalog#10438-H08H; NP_001314.1; Met1-Met149).
Ig Type:	Rabbit IgG
Clone ID:	022
Specificity:	Human Cystatin E / CST6
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA(Cap)
RECOMMENDED CONCENTRATION	
Sandwich ELISA (Capture Ab)	ELISA(Cap): 1:250-1:2000 This antibody will detect Human Cystatin E / CST6 in ELISA pair set (Catalog: # SEK10438). In a sandwich ELISA, it can be used as capture antibody when paired with (Catalog: # 10438-MM02).

Please Note: Optimal concentrations/dilutions should be determined by the end user.