

Carbonic Anhydrase XIV / CA14 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10458-R004

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human CA14 protein (Catalog#10458-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CA14 / Car14 extracellular domain (rh CA14; Catalog#10458-H08H; NP_036245.1; Met 1-Met 290).
Ig Type:	Rabbit IgG
Clone ID:	004
Specificity:	Human CA14 / Car14
	No cross-reactivity in ELISA with Human CAIX / CA9 Human CA10 / CARPX Human CA12 / CAXII Mouse CA14 / Car14
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human CA14.

Please Note: Optimal concentrations/dilutions should be determined by the end user.