

BID Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10468-T52

GENERAL INFORMATION

Immunogen:	Recombinant Human BID Protein (Catalog#10468-HNCE)
Preparation	Produced in rabbits immunized with purified, recombinant Human BID (rh BID; Catalog#10468-HNCE; P55957-1; Met1-Asp195). BID specific IgG was purified by Human BID affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human BID
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, ELISA, IP
---------------	---------------

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:2000
Immunoprecipitation	IP: 1-5 µL/mg of lysate
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

BID Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10468-T52




Anti-BID rabbit polyclonal antibody at 1:500 dilution
Lane A: A431 Whole Cell Lysate
Lane B: HepG2 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.
Secondary
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 22 kDa
Observed band size: 22 kDa



BID was immunoprecipitated using:
Lane A: 0.5 mg A431 Whole Cell Lysate
Lane B: 0.5 mg HepG2 Whole Cell Lysate

4 µL anti-BID rabbit polyclonal antibody and 60 µg of Immunomagnetic beads
Protein A/G.

Primary antibody:
Anti-BID rabbit polyclonal antibody, at 1:100 dilution

Secondary antibody:
Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 22 kDa
Observed band size: 22 kDa