

ATP5B Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 104727-T40

GENERAL INFORMATION	
Immunogen:	A synthetic peptide corresponding to the center region of the Human ATP5B
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human ATP5B, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human ATP5B
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB, IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:2500-1:10000
Western Blot	WB: 1:500-1:1000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

ATP5B Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 104727-T40




Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human ATP5B in human kidney with rabbit polyclonal antibody at 1:5000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human ATP5B in human duodenum with rabbit polyclonal antibody at 1:5000 dilution, formalin-fixed paraffin embedded sections.



Anti-ATP5B rabbit polyclonal antibody at 1:500 dilution

Lane A: HeLa Whole Cell Lysate

Lane B: HepG2 Whole Cell Lysate

Lane C: Mouse heart tissue lysate

Lane D: Rat heart tissue lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 56 kDa

Observed band size: 54 kDa

(We are unsure as to the identity of these extra bands.)