

CD131 / CSF2RB / IL3RB / IL5RB Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10516-MM11

GENERAL INFORMATION

Immunogen:	Recombinant Human CSF2RB protein (Catalog#10516-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD131 / CSF2RB / IL3RB / IL5RB (rh CD131 / CSF2RB / IL3RB / IL5RB; Catalog#10516-H08H; NP_000386.1; Met 1-Trp 443). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG2a
Clone ID:	5C7B1
Specificity:	Human CD131 / CSF2RB / IL3RB / IL5RB
	No cross-reactivity in ELISA with Human cell lysate (293 cell line)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human CSF2Rb.
--------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.