

CD208 / LAMP3 / DC-LAMP Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10527-RP01

GENERAL INFORMATION	
Immunogen:	Recombinant Human LAMP3 protein (Catalog#10527-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human CD208 / LAMP3 (rh LAMP3; Catalog#10527-H08H; NP_055213.2; Lys 28-Thr 381). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CD208 / LAMP3 / DC-LAMP
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human LAMP3.

Please Note: Optimal concentrations/dilutions should be determined by the end user.