

DOPA Decarboxylase / DDC Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10560-R005

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human DOPA Decarboxylase / DDC protein (Catalog#10560-H08B)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human DOPA Decarboxylase / DDC (rh DOPA Decarboxylase / DDC; Catalog#10560-H08B; NP_000781.1; Met 1-Glu 480).
Ig Type:	Rabbit IgG
Clone ID:	005
Specificity:	Human DOPA Decarboxylase / DDC
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA, IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:100-1:500
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human DDC.

Please Note: Optimal concentrations/dilutions should be determined by the end user.


DOPA Decarboxylase / DDC Antibody, Rabbit MAb

Catalog Number: 10560-R005



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



DDC Antibody	+
DDC Antigen	-
	+

Immunohistochemical staining of human DDC in human brain with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections). The left panel: tissue incubated with primary antibody; The right panel: tissue incubated with the mixture of primary antibody and antigen (recombinant protein).