

C5a Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10604-MM32

GENERAL INFORMATION

Immunogen:	10604-HNAE
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human C5a (rh C5a; Catalog#10604-HNAE; NP_001726.2; Leu679-Arg751). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	32
Specificity:	Human C5a
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	C5a

APPLICATIONS

Applications:	ELISA(Det)
---------------	------------

RECOMMENDED CONCENTRATION

Sandwich ELISA (Detection Ab)	ELISA(Det): 1:1000-1:10000 This antibody will detect Human C5a in ELISA pair set (Catalog: # SEKA10604). In a sandwich ELISA, it can be used as detection antibody when paired with (Catalog: # 10604-MM08).
-------------------------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>