

Complement C5a Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10604-R241

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	10604-HNAE
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human Complement C5a (rh Complement C5a ; Catalog#10604-HNAE; NP_001726.2; Leu679-Arg751).
Ig Type:	Rabbit IgG
Clone ID:	241
Specificity:	Human Complement C5a
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.
Alternative Names:	C5a
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000

Please Note: Optimal concentrations/dilutions should be determined by the end user.