

CAMKI Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 106746-T32

GENERAL INFORMATION	
Immunogen:	A synthetic peptide corresponding to the C-terminus of the Mouse CAMKI
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Mouse CAMKI, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Mouse Rat (Species predicted to react based on 100% sequence homology)
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:10000-1:50000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CAMKI Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 106746-T32



Anti-CAMK1 rabbit polyclonal antibody at
1:10000 dilution

Lane A: HepG2 Whole Cell Lysate

Lane B: Mouse brain tissue lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000
dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 41 kDa

Observed band size: 44 kDa