

# c-Met / HGFR Antibody (HRP), Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 10692-T16-H

## GENERAL INFORMATION

|                    |                                                                                                                                                                                                                                                                                                |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen:         | Recombinant Human c-Met / HGFR protein (Catalog#10692-H08H)                                                                                                                                                                                                                                    |
| Preparation        | Produced in rabbits immunized with purified, recombinant Human c-Met / HGFR (rh c-Met / HGFR; Catalog#10692-H08H; NP_000236.2; Met1-Thr932). c-Met / HGFR specific IgG was purified by Human c-Met / HGFR affinity chromatography and conjugated with horseradish-peroxidase (HRP).            |
| Ig Type:           | Rabbit IgG                                                                                                                                                                                                                                                                                     |
| Specificity:       | Human c-Met / HGFR                                                                                                                                                                                                                                                                             |
| Formulation:       | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                                                                |
| Storage:           | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles. |
| Alternative Names: | AUTS9                                                                                                                                                                                                                                                                                          |

## APPLICATIONS

|               |       |
|---------------|-------|
| Applications: | ELISA |
|---------------|-------|

## RECOMMENDED CONCENTRATION

|       |                    |
|-------|--------------------|
| ELISA | ELISA: 0.1-1 µg/mL |
|-------|--------------------|

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**

**For Research Use Only. Not for use in diagnostic or therapeutic procedures.**

**Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)**

**Website: <http://www.sinobiological.com>**