

EpCAM / TROP-1 / TACSTD1 Antibody (APC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10694-MM06-A

GENERAL INFORMATION

Immunogen:	Recombinant Human EpCAM / TROP-1 / TACSTD1 protein (Catalog#10694-H08H)
Reagents:	APC-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human EpCAM / TROP-1 / TACSTD1 (rh EpCAM / TROP-1 / TACSTD1; Catalog#10694-H08H; NP_002345.1; Met1-Lys265) and conjugated with APC under optimum conditions, the unreacted APC was removed.
Ig Type:	Mouse IgG1
Clone ID:	06
Specificity:	Human EpCAM / TROP-1 / TACSTD1
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !
Alternative Names:	EpCAM

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

EpCAM / TROP-1 / TACSTD1 Antibody (APC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10694-MM06-A



Profile of anti-EpCAM (CD326) reactivity on
SKBR3 cells analyzed by flow cytometry.

Flow cytometry was performed on a BD
FACSCalibur flow cytometry system.
Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for
technical protocols.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>