

CD86 / B7-2 / B70 Antibody (HRP), Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10699-R118-H

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human CD86 protein (Catalog#10699-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD86 / B7-2 (rh CD86; Catalog#10699-H08H; NP_008820.2; Met 1-His 239), and conjugated with horseradish-peroxidase (HRP).
Ig Type:	Rabbit IgG
Clone ID:	118
Specificity:	Human CD86 / B7-2 / B70
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Protected from prolonged exposure to light. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	CD86
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/ml

Please Note: Optimal concentrations/dilutions should be determined by the end user.