

IFN-beta Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10704-R865

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human IFN-beta Protein (Catalog#10704-HNAS)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human IFN-beta (rh IFN-beta; Catalog#10704-HNAS; P01574; Met1-Asn187).
Ig Type:	Rabbit IgG
Clone ID:	865
Specificity:	Human IFN-beta
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000
--------------	-----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.