



Catalog Number: 10707-R010

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human EPOR Protein (Catalog#10707-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human EPOR (rh EPOR; Catalog#10707-H08H; Met 1-Pro 250; NP_000112.1).
Ig Type:	Rabbit IgG
Clone ID:	010
Specificity:	Human EPOR
	Has cross-reactivity in ELISA with Mouse EPOR
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA,ELISA(Det)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human EPOR.
Sandwich ELISA (Detection Ab)	ELISA(Det): 1:1000-1:10000 This antibody will detect Human EPOR / Erythropoietin Receptor in ELISA pair set (Catalog: # SEK10707). In a sandwich ELISA, it can be used as detection antibody when paired with (Catalog: # 10707-R004).

Please Note: Optimal concentrations/dilutions should be determined by the end user.