

DPP7 / DPPII / DPP2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10748-RP04

GENERAL INFORMATION

Immunogen:	Recombinant Human DPPII protein (Catalog#10748-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human DPPII (rhDPPII; Catalog#10748-H08H; NP_037511.2, Met 1-Leu 492). DPPII specific IgG was purified by human DPPII affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human DPPII / DPP7
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human DPPII.
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.