

CD36 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10752-R331

EliteRmab® is a registered trademark of Sino Biological Inc.


GENERAL INFORMATION	
Immunogen:	Recombinant Human CD36 Protein (Catalog#10752-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CD36 (rh CD36; Catalog#10752-H08H; NP_001001547.1; Gly30-Asn439).
Ig Type:	Rabbit IgG
Clone ID:	331
Specificity:	Human CD36
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	BDPLT10,CHDS7,FAT,GP3B,GP4,GPIV,PASIV,SCARB3
APPLICATIONS	
Applications:	IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:100-1:500

Please Note: Optimal concentrations/dilutions should be determined by the end user.




Catalog Number: 10752-R331


EliteRmab® is a registered trademark of Sino Biological Inc.



Immunohistochemical staining of human CD36 in human liver with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CD36 in human stomach with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human CD36 in human placenta with rabbit monoclonal antibody at 1:200 dilution, formalin-fixed paraffin embedded sections.