

Adiponectin Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10772-MM12

GENERAL INFORMATION

Immunogen:	Recombinant Human Adiponectin protein (Catalog#10772-H08B)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human Adiponectin (rh Adiponectin; Catalog#10772-H08B; NP_001171271.1; Met1-Asn244). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Clone ID:	12
Ig Type:	mouse IgG1
Specificity:	Human Adiponectin
Formulation:	0.2 µm filtered solution in PBS with 10% Trehalose, pH7.0
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA,ELISA(Det)
---------------	------------------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000
ELISA(Det)	ELISA(Det): 1:1000-1:10000

Please Note: Optimal concentrations/dilutions should be determined by the end user.