

ASGPR1 / ASGR1 Antibody (FITC), Mouse MAbs



Sino Biological
Biological Solution Specialist

Catalog Number: 10773-MM02-F

GENERAL INFORMATION

Immunogen:	Recombinant Human ASGPR1 / ASGR1 Protein (Catalog#10773-H07H)
Reagents:	FITC-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human ASGPR1 / ASGR1 (rh ASGPR1 / ASGR1; Catalog#10773-H07H; P07306; Gln62-Leu291) and conjugated with FITC under optimum conditions, the unreacted FITC was removed.
Ig Type:	Mouse IgG1
Clone ID:	02
Specificity:	Human ASGPR1 / ASGR1
Concentration:	10 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !
Alternative Names:	ASGR1

APPLICATIONS

Applications:	FCM
---------------	-----

RECOMMENDED CONCENTRATION


Please Note: Optimal concentrations/dilutions should be determined by the end user.

ASGPR1 / ASGR1 Antibody (FITC), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10773-MM02-F



Flow cytometric analysis of Human Asialoglycoprotein Receptor expression on HepG2 cells. Cells were stained with FITC-conjugated anti-Human Asialoglycoprotein Receptor. The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.