

Asialoglycoprotein Receptor Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10773-R009

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human Asialoglycoprotein Receptor Protein (Catalog#10773-H07H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human Asialoglycoprotein Receptor (rh Asialoglycoprotein Receptor; Catalog#10773-H07H; P07306; Gln62-Leu291).
Ig Type:	Rabbit IgG
Clone ID:	009
Specificity:	Human Asialoglycoprotein Receptor
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ICC/IF
----------------------	--------

RECOMMENDED CONCENTRATION

ICC/IF	ICC/IF: 1:20-1:100
---------------	--------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.


Asialoglycoprotein Receptor Antibody, Rabbit MAb

Catalog Number: 10773-R009



Sino Biological
Biological Solution Specialist

EliteRmab® is a registered trademark of Sino Biological Inc.



Immunofluorescence staining of ASGR1 in HepG2 cells. Cells were fixed with 4% PFA, permeabilized with 0.1% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-Human ASGR1 monoclonal antibody (dilution ratio 1:60) at 4°C overnight. Then cells were stained with the Alexa Fluor®488-conjugated Goat Anti-rabbit IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to Cytoplasm.