



Catalog Number: 10773-R011

EliteRmab® is a registered trademark of Sino Biological Inc.


GENERAL INFORMATION	
Immunogen:	Recombinant Human ASGPR1 protein (Catalog#10773-H07H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human ASGPR1 (rh ASGPR1 ; Catalog#10773-H07H; P07306; Gln62-Leu291).
Ig Type:	Rabbit IgG
Clone ID:	011
Specificity:	Human ASGPR1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	IHC-P
RECOMMENDED CONCENTRATION	
IHC-P	IHC-P: 1:100-1:500

Please Note: Optimal concentrations/dilutions should be determined by the end user.




Catalog Number: 10773-R011

EliteRmab® is a registered trademark of Sino Biological Inc.



Immunochemical staining of human ASGR1 in human hepatoma with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).



Immunochemical staining of human ASGR1 in human cirrhosis with rabbit monoclonal antibody (1:200, formalin-fixed paraffin embedded sections).