

CD48 / Blast-1 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10797-MM04

GENERAL INFORMATION

Immunogen:	Recombinant Human CD48 / Blast-1 Protein (Catalog#10797-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CD48 / Blast-1 (rh CD48 / Blast-1; Catalog#10797-H08H; NP_001769.2; Met1-Ser220). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Clone ID:	04
Ig Type:	Mouse IgG1
Specificity:	Human CD48 / Blast-1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
---------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000
-------	----------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.