

GSS/Glutathione Synthetase Antibody, Rabbit PAb, Antigen Affinity Purified

Catalog Number: 108030-T08



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the center region of the Human GSS/Glutathione Synthetase
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the center region of the Human GSS/Glutathione Synthetase, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human GSS/Glutathione Synthetase
Formulation:	PBS, pH7.0 with 0.03% Proclin300
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
----------------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:500-1:2000
--------------	---------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

**GSS/Glutathione Synthetase
Antibody, Rabbit PAb, Antigen
Affinity Purified**


Catalog Number: 108030-T08



Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human GSS in human kidney with rabbit polyclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human GSS in human liver with rabbit polyclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.