

CD38 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10818-T58

GENERAL INFORMATION

Immunogen:	Recombinant Human CD38 Protein (Catalog#10818-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human CD38 (rh CD38; Catalog#10818-H08H; NP_001766.2; Val43-Ile300). CD38 specific IgG was purified by Human CD38 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human CD38
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, ELISA, IHC-P, ICC/IF
----------------------	--------------------------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:500-1:2000
Western Blot	WB: 1:500-1:2000
IF	ICC/IF: 1:50-1:200
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

CD38 Antibody, Rabbit PAb, Antigen Affinity Purified


Catalog Number: 10818-T58




Sino Biological
Biological Solution Specialist



Immunohistochemical staining of human CD38 in human tonsil with rabbit polyclonal antibody at 1:1000 dilution, formalin-fixed paraffin embedded sections.



Immunofluorescence staining of CD38-His in Raji cells. Cells were fixed with 4% PFA, blocked with 10% serum, and incubated with rabbit anti-Human CD38-His polyclonal antibody (dilution ratio 1:100) at 4°C overnight. Then cells were stained with the Alexa Fluor®488-conjugated Goat Anti-rabbit IgG secondary antibody (green) and counterstained with DAPI (blue). Positive staining was localized to Cell membrane.



Anti-CD38 rabbit polyclonal antibody at 1:500 dilution

Lane A: Jurkat Whole Cell Lysate

Lane B: HuT-78 Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.
Performed under reducing conditions.

Predicted band size: 34 kDa

Observed band size: 38 kDa