

CER1 / Cerberus 1 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10820-R034

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human CER1 / Cerberus 1 Protein (Catalog#10820-H08B)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human CER1 / Cerberus 1 (rh CER1 / Cerberus 1; Catalog#10820-H08B; AAH69371.1; Met1-Ala267).
Ig Type:	Rabbit IgG
Clone ID:	034
Specificity:	Human CER1 / Cerberus 1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human CER1 / Cerberus 1.

Please Note: Optimal concentrations/dilutions should be determined by the end user.