

CEACAM6 / CD66c Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10823-MM12-P

GENERAL INFORMATION

Immunogen:	Recombinant Human CEACAM6 / CD66c protein (Catalog#10823-H08H)
Reagents:	PE-conjugated Mouse monoclonal antibody
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CEACAM6 / CD66c (rh CEACAM6 / CD66c; Catalog#10823-H08H; NP_002474.4; Met1-Gly320) and conjugated with PE under optimum conditions, the unreacted PE was removed.
Ig Type:	Mouse IgG1
Clone ID:	12
Specificity:	Human CEACAM6 / CD66c
Concentration:	5 µl/Test, 0.1 mg/ml
Formulation:	PBS solution containing 0.5% BSA and 0.03%ProClin300
Storage:	This antibody can be stored at 2°C-8°C for twelve months without detectable loss of activity. Protected from prolonged exposure to light. Do not freeze !
Alternative Names:	CEACAM6

APPLICATIONS

Applications:	FCM
----------------------	-----

RECOMMENDED CONCENTRATION

Please Note: Optimal concentrations/dilutions should be determined by the end user.

CEACAM6 / CD66c Antibody (PE), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10823-MM12-P

Flow cytometric analysis of human CEACAM6 expression on human whole blood granulocytes. Human whole blood granulocytes were stained with PE-conjugated anti-Human CEACAM6. The histogram were derived from gated events with the forward and side light-scatter characteristics of viable granulocytes.

Flow cytometry was performed on a BD FACSCalibur flow cytometry system. Please refer to www.sinobiological.com/Flow-Cytometry-FACS-Protocols-a-750.html for technical protocols.