

FGFR2 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10824-MM02T

GENERAL INFORMATION	
Immunogen:	Recombinant Human FGFR2 Protein (Catalog#10824-H08H)
Preparation	This product is a recombinant monoclonal antibody expressed from HEK293 cells.
Clone ID:	02
Specificity:	Human FGFR2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
Alternative Names:	BBDS, BEK, BFR-1, CD332, CEK3, CFD1, ECT1, JWS, K-SAM, KGFR, TK14, TK25
APPLICATIONS	
Applications:	ELISA(Det)
RECOMMENDED CONCENTRATION	
Sandwich ELISA (Detection Ab)	ELISA(Det): 1:1000-1:10000 This antibody will detect Human FGFR2 / CD332 in ELISA pair set (Catalog: # SEK10824). In a sandwich ELISA, it can be used as detection antibody when paired with (Catalog: # 10824-MM10T).

Please Note: Optimal concentrations/dilutions should be determined by the end user.