

FGFR2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10824-T48

GENERAL INFORMATION

Immunogen:	Recombinant Human FGFR2 Protein (Catalog#10824-H08H)
Preparation	Produced in rabbits immunized with purified, recombinant Human FGFR2 (rh FGFR2; Catalog#10824-H08H; NP_000132.3; Met1-Glu377). FGFR2 specific IgG was purified by Human FGFR2 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human FGFR2
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB, ELISA
---------------	-----------

RECOMMENDED CONCENTRATION

Western Blot	WB: 1:500-1:2000
ELISA	ELISA: 1:5000-1:10000


Please Note: Optimal concentrations/dilutions should be determined by the end user.

FGFR2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10824-T48



Anti-FGFR2 rabbit polyclonal antibody at 1:500 dilution
Lane A: Jurkat Whole Cell Lysate

Lysates/proteins at 30 µg per lane.

Secondary

Goat Anti-Rabbit IgG (H+L)/HRP at 1/10000 dilution.

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 92 kDa

Observed band size: 150 kDa

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Tel: +86-400-890-9989 (Global), +1-215-583-7898 (USA), +49(0)6196 9678656 (Europe)

Website: <http://www.sinobiological.com>