

Hemopexin / HPX Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10870-MM02

GENERAL INFORMATION	
Immunogen:	Recombinant Human HPX protein (Catalog#10870-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human HPX / Hemopexin (rhHPX; Catalog#10870-H08H; Met 1-His 462; NP_000604.1). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	2H4C4F12
Specificity:	Human HPX / Hemopexin
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human HPX.

Please Note: Optimal concentrations/dilutions should be determined by the end user.