

I-TAC / CXCL11 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10876-MM12

GENERAL INFORMATION	
Immunogen:	Recombinant Human I-TAC / CXCL11 Protein (Catalog#10876-HNAE)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human I-TAC / CXCL11 (rh I-TAC / CXCL11; Catalog#10876-HNAE; O14625; Phe22-Phe94). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG2a
Clone ID:	12
Specificity:	Human I-TAC / CXCL11
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human I-TAC / CXCL11.

Please Note: Optimal concentrations/dilutions should be determined by the end user.