

Cytokeratin 8 / KRT8 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 108772-T08

GENERAL INFORMATION

Immunogen:	A synthetic peptide corresponding to the C-terminus of the Human Cytokeratin 8 / KRT8
Preparation	Produced in rabbits immunized with a synthetic peptide corresponding to the C-terminus of the Human Cytokeratin 8 / KRT8, and purified by antigen affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human Cytokeratin 8 / KRT8
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	IHC-P
---------------	-------

RECOMMENDED CONCENTRATION

IHC-P	IHC-P: 1:200-1:1000
-------	---------------------


Please Note: Optimal concentrations/dilutions should be determined by the end user.

Cytokeratin 8 / KRT8 Antibody, Rabbit PAb, Antigen Affinity Purified




Sino Biological
Biological Solution Specialist


Catalog Number: 108772-T08



Immunohistochemical staining of human KRT8 in human colon carcinoma with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human KRT8 in human prostate with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.



Immunohistochemical staining of human KRT8 in human small intestine with rabbit polyclonal antibody at 1:300 dilution, formalin-fixed paraffin embedded sections.