

Ephrin-A1 / EFNA1 Antibody, Rabbit MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10882-R114

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
Immunogen:	Recombinant Human EFNA1 protein (Catalog#10882-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human EphrinA1 (rh EphrinA1; Catalog#10882-H08H; NP_004419.2; Met 1-Ser 182).
Ig Type:	Rabbit IgG
Clone ID:	114
Specificity:	Human EphrinA1 / EFNA1
	No cross-reactivity in ELISA with Human EphrinA5 / EFNA5 Human EphrinA3 / EFNA3
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA,ELISA(Det)
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human EFNA1.
Sandwich ELISA (Detection Ab)	ELISA(Det): 1:1000-1:10000 This antibody will detect Human Ephrin-A1 / EFNA1 in ELISA pair set (Catalog: # SEK10882). In a sandwich ELISA, it can be used as detection antibody when paired with (Catalog: # 10882-MM01).

Please Note: Optimal concentrations/dilutions should be determined by the end user.