

CXCL5 Antibody (Biotin), Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10889-MM14-B

GENERAL INFORMATION	
Immunogen:	Recombinant Human CXCL5 Protein (Catalog#10889-HNAE)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human CXCL5 (rh CXCL5; Catalog#10889-HNAE; P42830; Arg45-Asn114). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography and then biotinylated.
Ig Type:	Mouse IgG1
Clone ID:	14
Specificity:	Human CXCL5
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	ELISA
RECOMMENDED CONCENTRATION	
ELISA	ELISA: 0.1-1 µg/mL

Please Note: Optimal concentrations/dilutions should be determined by the end user.